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




<p>1. </p> <p>LABEL microcentrifuge tubes A-E for the five CRISPR-Cas9 digestions. Put your initials or group number on the tube.</p>	<p>2. </p> <p>Use a micropipette to ADD 15 μL of SH2 DNA to tubes A-E.</p>	<p>3. </p> <p>ADD 10 μL of Cas9 to tubes A-E.</p>	<p>4. </p> <p>Use Table 2 below to ADD 15 μL of each gRNA to its respective tube.</p>
<p>5. </p> <p>CAP tubes A-E and TAP gently on the lab bench to mix and collect contents at the bottom of the tubes.</p>	<p>6. </p> <p>INCUBATE tubes in a 37°C water bath for 30-60 min. (Alternatively, tubes can be incubated in a 45°C water bath for 15 min.).</p>	<p>7. </p> <p>ADD 5 μL of Gel Loading Solution to tubes A-E to stop the reactions.</p>	<p>8. </p> <p>CAP tubes A-E and TAP gently on the lab bench to mix. PROCEED to Module IV: Agarose Gel Electrophoresis.</p>

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<p>1. </p> <p>DILUTE concentrated (50X) electrophoresis buffer with distilled water to create 1X buffer (see Table A).</p>	<p>2. </p> <p>MIX agarose powder with 1X buffer in a flask (see Table A).</p>	<p>3. </p> <p>DISSOLVE agarose powder by boiling the solution. MICROWAVE the solution on high for 1 minute. Carefully REMOVE the flask and MIX by swirling. Continue to HEAT the solution in 15-second bursts until the agarose is completely dissolved (the solution should be clear like water).</p>
<p>4. </p> <p>While agarose is cooling, SEAL the ends of the gel-casting tray with the rubber end caps. PLACE the well template (comb) in the appropriate notch.</p>	<p>5. </p> <p>Before casting the gel, ADD diluted SYBR® Safe to the molten agarose (see Table A) and SWIRL to mix.</p>	<p>6. </p> <p>POUR the cooled agarose solution into the prepared gel-casting tray. The gel should thoroughly solidify within 20 minutes. The gel will stiffen and become less transparent as it solidifies.</p>

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


Size of Gel Casting Tray	Concentrated Buffer (50x)	+ Distilled Water	+ Amt of Agarose	= TOTAL Volume	Diluted SYBR® (Step 6)
7 x 7 cm	0.6 mL	29.4 mL	0.24 g	30 mL	30 µL
10 x 7 cm*	0.9 mL	44.1 mL	0.36 g	45 mL	45 µL
14 x 7 cm	1.2 mL	58.8 mL	0.48 g	60 mL	60 µL

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<p>8.</p>  <p>REMOVE end caps and comb. Take particular care when removing the comb to prevent damage to the wells.</p>	<p>9.</p>  <p>PLACE the gel (still on the tray) into the electrophoresis chamber. COVER the gel with 1X Electrophoresis Buffer (See Table B for recommended volumes). The gel should be completely submerged.</p>	<p>10.</p>  <p>OBTAIN your sample tubes. Make sure the gel is properly oriented in the chamber. LOAD the entire sample (35 μL) into the well in the order indicated by Table 3, below.</p> <table border="1" style="margin-left: auto; margin-right: auto; border-collapse: collapse; text-align: center;"> <thead> <tr> <th colspan="2">TABLE 3: GEL LOADING</th> </tr> </thead> <tbody> <tr> <td>Lane 1</td> <td>DNA Standard Marker</td> </tr> <tr> <td>2</td> <td>Tube A - gRNA A</td> </tr> <tr> <td>3</td> <td>Tube B - gRNA B</td> </tr> <tr> <td>4</td> <td>Tube C - gRNA C</td> </tr> <tr> <td>5</td> <td>Tube D - gRNA D</td> </tr> <tr> <td>6</td> <td>Tube E - Scramble gRNA</td> </tr> </tbody> </table>	TABLE 3: GEL LOADING		Lane 1	DNA Standard Marker	2	Tube A - gRNA A	3	Tube B - gRNA B	4	Tube C - gRNA C	5	Tube D - gRNA D	6	Tube E - Scramble gRNA
TABLE 3: GEL LOADING																
Lane 1	DNA Standard Marker															
2	Tube A - gRNA A															
3	Tube B - gRNA B															
4	Tube C - gRNA C															
5	Tube D - gRNA D															
6	Tube E - Scramble gRNA															
<p>11.</p>  <p>PLACE safety cover on the unit. CHECK that the gel is properly oriented. Remember, the DNA samples will migrate toward the positive (red) electrode.</p>	<p>12.</p> <p>CONNECT leads to the power source and PERFORM electrophoresis (See Table C for time and voltage guidelines). Allow the tracking dye to migrate at least 3.5 cm from the wells.</p> 	<p>13.</p> <p>After electrophoresis is complete, REMOVE the gel and casting tray from the electrophoresis chamber and proceed to Step 14 for VISUALIZING the agarose gel.</p>														

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



EDVOTEK Model #	Total Volume Required	Dilution 50x Conc. Buffer + Distilled Water	
EDGE™	150 mL	3 mL	147 mL
M12	400 mL	8 mL	392 mL
M36	1000 mL	20 mL	980 mL

Volts	Electrophoresis Model	
	EDGE™	M12 & M36
	Min/Max (minutes)	Min/Max (minutes)
150	10/20	20/35
125	N/A	30/45
100	15/25	40/60

<p>14.</p>  <p>SLIDE gel off the casting tray onto the viewing surface of the transilluminator. LOWER orange contrast cover.</p>  <p>Be sure to wear UV goggles if using a UV transilluminator.</p>	<p>15.</p>  <p>Turn the unit ON. DNA should appear as bright green bands on a dark background.</p> <p>PHOTOGRAPH the results.</p>	<p>16.</p>  <p>Turn the unit OFF. REMOVE and DISPOSE of the gel. CLEAN the transilluminator surfaces with distilled water.</p>
---	---	---

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
Sample Name	Top Hit from BLAST
gRNA A	<i>Zea mays</i> ADP-glucose pyrophosphorylase large subunit (shrunken2) - SH2
gRNA B	<i>Zea mays</i> shrunken-2 (SH2) gene locus, partial sequence
gRNA C	<i>Zea mays</i> shrunken-2 (SH2) gene locus, partial sequence
gRNA D	<i>Saccharum hybrid</i> cultivar R570 done BAC 022D05 complete sequence

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MODULE IV

Lane	Sample	Result	Molecular Weights
1	Standard DNA Marker	---	6751, 3652, 2877, 1568, 1118, 825, 630
2	Tube A - gRNA A	DNA is cut	3000, 1280
3	Tube B - gRNA B	Uncut DNA	4280
4	Tube C - gRNA C	DNA is cut	3000, 1280
5	Tube D - gRNA D	DNA is cut	3650, 630
6	Tube E - Scramble gRNA	Uncut DNA	4280