Vernier Photogate 포토게이트센서

Order Code VPG-BTD

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Blocked To Blocked	<p>("Time", "Gate1", "Gate2") "Time": Optional. A column of real numbers (the times of events) "Gate1": A column of photogate states (1's and 0's) "Gate2": Optional. A column of photogate states (1's and 0's) Returns a column of the times between successive blocked events in gate 1 and blocked events in gate 2. If you don't enter a "Time" column, the program will find one. If you don't enter "Gate2", "Gate1" will be used.</p>
Blocked to Unblocked	<p>("Time", "Gate1", "Gate2") "Time": Optional. A column of real numbers (the times of events) "Gate1": A column of photogate states (1's and 0's) "Gate2": Optional. A column of photogate states (1's and 0's) Returns a column of the times between successive blocked events in gate 1 and unblocked events in gate 2. If you don't enter a "Time" column, the program will find one. If you don't enter "Gate2", "Gate1" will be used.</p>
Unblocked to Blocked	<p>("Time", "Gate1", "Gate2") "Time": Optional. A column of real numbers (the times of events) "Gate1": A column of photogate states (1's and 0's) "Gate2": Optional. A column of photogate states (1's and 0's) Returns a column of the times between successive unblocked events in gate 1 and blocked events in gate 2. If you don't enter a "Time" column, the program will find one. If you don't enter "Gate2", "Gate1" will be used.</p>
Unblocked to Unblocked	<p>("Time", "Gate1", "Gate2") "Time": Optional. A column of real numbers (the times of events) "Gate1": A column of photogate states (1's and 0's) "Gate2": Optional. A column of photogate states (1's and 0's) Returns a column of the times between successive unblocked events in gate 1 and unblocked events in gate 2. If you don't enter a "Time" column, the program will find one. If you don't enter "Gate2", "Gate1" will be used.</p>
Blocked Mid Times	<p>("Time", "Gate1", "Gate2") "Time": Optional. A column of real numbers (the times of events) "Gate1": A column of photogate states (1's and 0's) "Gate2": Optional. A column of photogate states (1's and 0's) Calculate the average times between blocked events from Gate 1 to Gate 2. If you don't enter a "Time" column, the program will find one. If you don't enter "Gate2", "Gate1" will be used.</p>
Blocked to Unblocked MidTimes	<p>Blocked to Unblocked MidTimes "Time": Optional. A column of real numbers (the times of events) "Gate1": A column of photogate states (1's and 0's) "Gate2": Optional. A column of photogate states (1's and 0's) Calculate the average time between successive blocked events in gate 1 and unblocked events in gate 2. If you don't enter a "Time" column, the program will find one. If you don't enter "Gate2", "Gate1" will be used.</p>
Unblocked to Blocked MidTimes	<p>Unblocked to Blocked MidTimes "Time": Optional. A column of real numbers (the times of events) "Gate1": A column of photogate states (1's and 0's) "Gate2": Optional. A column of photogate states (1's and 0's) Calculate the average time between successive unblocked events in gate 1 and blocked events in gate 2. If you don't enter a "Time" column, the program will find one. If you don't enter "Gate2", "Gate1" will be used.</p>

Unblocked Mid Times	<p>("Time", "Gate1", "Gate2") "Time": Optional. A column of real numbers (the times of events) "Gate1": A column of photogate states (1's and 0's) "Gate2": Optional. A column of photogate states (1's and 0's) Calculate the average times between blocked events from Gate 1 to Gate 2. If you don't enter a "Time" column, the program will find one. If you don't enter "Gate2", "Gate1" will be used.</p>
Pendulum Period	<p>("Time", "Gate1") "Time": Optional. A column of real numbers (the times of events) "Gate1": A column of photogate states (1's and 0's) Calculate the time between every other blocked event on Gate 1. If you don't enter a "Time" column, the program will find one.</p>
derivative Time Shift	<p>DerivativeTimeShift ("Y", "X"): returns the derivative of "Y" with respect to "X". This function is specifically designed to be used with photogate and picket fence data. The derivatives returned are adjusted to estimate values at the start of the timing interval, instead of the midpoint. For details see The Physics Teacher, Vol 35, April 1997, p.220. Average velocity during the time interval is equal to the instantaneous velocity a $v_2 = \frac{v_{ave,12} \Delta t_{12} + v_{ave,23} \Delta t_{12}}{\Delta t_{13}}$ val.</p> <p>Where $v_{ave,12} = \Delta x_{12} / \Delta t_{12}$ and $v_{ave,23} = \Delta x_{23} / \Delta t_{23}$</p>
second Derivative Time Shift	<p>Numerical time-shifted second derivative("Y", "X") "Y": A column of real numbers "X": Optional. A column of real numbers Numerical time-shifted second derivative. Calculates the second numerical derivative of "Y" with respect to "X". The values are shifted so that the derivative $a_2 = 2 \frac{v_{ave,23} - v_{ave,12}}{t_3 - t_1}$ points between each two values. If you don't supply a will find one.</p>
collapse	<p>collapse("X"): Returns a column with all the non-numerical cells (blanks and text) removed.</p>
collapse Indirect	<p>collapse Indirect(X, Y) Returns a column of only the rows in "X" corresponding to rows in "Y" that have valid numerical cells.</p>

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